KIT EVEREST

855.SET

Regolamento (UE)
2017/745

IT

WWW.KONG.IT

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CAPITOLO 9

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9.1 ETICHETTATURA DEL DM

Label for DM 855.00 - MOD. EVEREST. The label contains the KONG logo, a QR code, and various identification codes: MD, REF, UDI, and SN (123456 1234). It also features the CE mark, an information icon, and a factory icon. The text includes the model name, a UDI number (01) 08023577044446(21) 1234561234, and the manufacturer's address: KONG S.p.A., Via XXV Aprile, 4, I - 23804 MONTE MARENZO (LC) - ITALY, Tel. +39 0341 630506. The origin is marked as ITALY 2022 / 01.

Label for DM 855.11 - MOD. X-TRIM 2. The label contains the KONG logo, a QR code, and various identification codes: MD, REF, UDI, and SN (123456 1234). It also features the CE mark, an information icon, and a factory icon. The text includes the model name, a UDI number (01) 08023577015057(21) 1234561234, and the manufacturer's address: KONG S.p.A., Via XXV Aprile, 4, I - 23804 MONTE MARENZO (LC) - ITALY, Tel. +39 0341 630506. The origin is marked as ITALY 2022 / 01.

Label for DM 855.06 - MOD. HRP. The label contains the KONG logo, a QR code, and various identification codes: MD, REF, UDI, and SN (123456 1234). It also features the CE mark, an information icon, and a factory icon. The text includes the model name, a UDI number (01) 08023577044606(21) 1234561234, and the manufacturer's address: KONG S.p.A., Via XXV Aprile, 4, I - 23804 MONTE MARENZO (LC) - ITALY, Tel. +39 0341 630506. The origin is marked as ITALY 2022 / 01.

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TABLE OF CONTENTS

1 - SYMBOLS AND SUPPORT	4
1.1 Symbols	4
1.2 Technical assistance	4
2 - GENERAL INFORMATION	5
3 - TECHNICAL FEATURES	6
3.1 Part terminology and materials	6
3.2 Dimensions	7
3.3 Capacity	7
3.4 Optional components and spare parts	8
4 - SPECIFIC INFORMATION	10
4.1 Application	10
4.2 Assembly of the spinal board	10
4.3 Mounting the head immobiliser on the spinal board	11
4.4 Immobilising the patient on the spinal board	12
4.5 Unlocking the patient	15
4.6 Positioning the spinal board inside the rescue bag	16
4.7 Transporting the patient	22
5 - MAINTENANCE AND REPAIR	23
5.1 General Information	23
5.2 Maintenance	23
5.3 Repair	23
6 - STORAGE	24
7 - CHECKS, INSPECTIONS AND REVISION	25
7.1 Checks	25
7.2 Inspections	25
7.3 Revisions	25
8 - DURATION AND WARRANTY	26
8.1 Duration	26
8.2 Disposal	26
8.3 Warranty	26
8.4 Law obligations	26
9 - LABELLING AND SYMBOLS	27
9.1 Labelling	27
9.1 Separable components label	28
9.2 Symbols	28
10 - DOCUMENTS	29
10.1 Inspections and servicing register	29
10.2 Maintenance and repair register	32
10.3 EU Declaration of conformity	35
11 - LIST OF APPLIED AND REFERENCE STANDARDS	38
11.1 Applied standards	38
11.2 Standards used as reference	39

CHAPTER

1

SYMBOLS AND SUPPORT

1.1 SYMBOLS

For the sake of comfort and clarity while reading this manual, below are the symbols used for handling important warnings for a proper and safe use of the device.

REQUIREMENT FOR PROPER USE

It identifies the presence of information for proper use of the device.

INFORMATION REQUIREMENT

It identifies the presence of useful and general information which guides the reader towards a conscious use of the device and/or performance of actions.

It identifies that the Medical Device is manufactured, designed and produced in compliance with the provisions of the General Safety and Performance Requirements of Regulation (EU) 2017/745 (Class I medical device, in compliance with classification rule 1 as indicated in the attachment VIII).

1.2 SUPPORT

For information please contact Kong Customer Support Service by:

- telephone +39 0341 630506
- fax +39 0341 641550
- email: safetycare@kong.it

or write to KONG S.p.A. – Via XXV Aprile, 4 – 23804 Monte Marenzo LC - ITALY.

To facilitate support operations, please always communicate or state the serial number (SN) indicated on the label fixed to the Medical Device.

CHAPTER 2

2 GENERAL INFORMATION

The information provided by the manufacturer (hereinafter referred to as information) must be read and clearly understood by the user prior to using the Medical Device (MD hereafter). The information regards the description of the features, performance, assembly, disassembly, maintenance, preservation, disinfection, etc. of the device. Even though the information offers tips on use, this information shall not be considered as a user manual under actual conditions of use.

WARNINGS AND LIMITATIONS OF USE:

- this device shall be strictly used by people who are physically fit, trained (instructed and taught) to use the device and with specific experience regarding moving the patient or, during training activities, by people under direct supervision of the trainers/ supervisors who guarantee the safety thereof,
- do not use the device before fully reading and understanding this user manual,
- strictly follow the manufacturer's information, improper use of the device is hazardous,
- modifying and/or repairing the device is strictly forbidden,
- all checks described in chapter 7 must be carried out prior to and after using the device. In case of any doubt on the efficiency of the device, the user must replace it immediately,
- non-compliant use, deformations, falls, wear, chemical contamination, exposing textile/plastic components/devices to temperature below -30°C or above

- +50°C and metal components/devices to temperatures exceeding 100°C are some examples of causes that can reduce, limit and end the life of the device,
- Prior to any rescue operation, be keen not to exceed the capacity indicated in paragraph 3.3,
- in order to reduce risks of exposure to / transmission of infectious diseases, clean and disinfect the device as indicated in chapter 5,
- improper use of the patient immobilisation systems can jeopardise the safety of the patient,
- always check the compatibility of the devices used alongside the device by consulting the relative manufacturer's information,
- use of spare parts or optional components different from the ones indicated in paragraph 3.4 can be hazardous,
- do not expose the device to sources of heat and at contact with chemical substances. Reduce exposure to direct sunlight as much as possible. At low temperatures and in the presence of humidity, formation of ice could reduce flexibility and increase cutting and abrasion-related risks on textile and synthetic devices,
- Report to the manufacturer and the competent authority of the member state where the user and / or patient is established any serious incident occurring in relation to the MD.

All our devices are tested/ checked piece by piece in compliance with the procedures laid down by the Quality System certified in accordance with the UNI EN ISO 9001 standard. Laboratory tests, testing, information and standards do not always reproduce the practical result. Thus, the results obtained under the actual conditions of use of the device in the natural environment may differ, even considerably at times. The best information lies in the continuous practical use under the supervision of skilled/ expert/qualified people.

Warning: not suitable for use in an ATEX environment (Directive 94/9/CE)

CHAPTER 3

TECHNICAL FEATURES

3.1 TERMINOLOGY AND MATERIALS OF THE PARTS

- A – Everest bag
- B – X-TRIM 2 Spinal board
- C – Taylan Belt
- D – Taylan X-Tension webbing
- E – Patty hanging kit
- F – HRP Head immobiliser
- G – Carriage backpack


Fig.1

3.2 DIMENSIONS

Overall length: 182 cm

Overall width: 42 cm

Overall weight: 14,5 kg

3.3 CAPACITY

The load capacity of the “EVEREST KIT” is defined according to the results of the laboratory tests and the recommended safety coefficients.

How to use	Capacity
Manual transportation (fig. 2)	150 kg
Lifting and lowering using rescue devices (fig. 3) <i>with 855.03 PATTY hanging kit</i>	
Lifting and lowering using helicopters (fig. 3) <i>with 855.03 PATTY hanging kit</i>	

Prior to any rescue operation, ensure not to exceed the indicated capacity.


Fig.2


Fig.3

3.4 OPTIONAL COMPONENTS AND SPARE PARTS

3.4.1 Optional components

855.09 EXPRESS FRAME

Light and modular frame that can be disassembled to easily move the Everest kit on shoulders.


Fig.4

871.14 VAR

Manual control anti-rotation system designed for hoist recovery of the Kit EVEREST.


Fig.5

3.4.2 Spare parts

- A – Everest bag
- B – X-TRIM 2 Spinal board
- C – Taylan Belt
- D – Taylan X-Tension webbing
- E – Patty hanging kit
- F – HRP Head immobiliser
- G – Carriage backpack


Fig.6

CHAPTER 4

SPECIFIC INFORMATION

4.1 APPLICATION

The “EVEREST KIT” is a device particularly suitable for helicopter rescue. Decisions concerning moving and immobilising the patient, as well as the duration, method to be used and combination with other devices, shall be taken and executed by expert and trained personnel only.

4.2 ASSEMBLY OF THE SPINAL BOARD

a) Open the spinal board (fig. 7)


Fig.7

b) Slide the stiffening spars until they lock into the lateral grooves of the spinal board (fig. 8)


Fig.8

4.3 MOUNTING THE HEAD IMMOBILISER ON THE SPINAL BOARD

a) Insert the Velcro straps into the slots of the head immobiliser and spinal board (fig. 9).


Fig.9

b) Attach the head immobiliser by pulling and overlapping the Velcro straps on the underside of the spinal board (fig. 10).


Fig.10

4.4 IMMOBILISING THE PATIENT ON THE SPINAL BOARD

With the spinal board positioned on a flat and stable surface:

- a) place the patient's head in the head immobiliser and immobilise it using the padded straps (figs. 11 e 12),

The head immobiliser is also suitable for immobilising patients with cervical collars.


Fig.11

Fig.12

- b) Position the “TAYLAN” belt as shown in fig. 13,


Fig.13

c) fasten the red straps (shoulders and torso) by hooking the connectors to the slots on the board, as shown in figs. 14 e 15, and tighten them slightly by using the buckles (figs. 16 e 17),


Fig.14


Fig.15


Fig.16


Fig.17

d) fasten the black straps (pelvis and legs) by hooking the connectors to the slots on the board and slightly tighten them using the buckles (figs. 18 e 19)


Fig.18


Fig.19

e) conclude the immobilisation by using the longitudinal lock (fig. 20),


Fig.20

f) the patient fully immobilised on the spinal board can be moved by the operators. Avoid prolonged contact of the patient's skin with the device's materials.

4.5 UNLOCKING THE PATIENT

a) Loosen the straps by using the buckles as shown in fig. 21.

If the straps are removed from their buckles, they can be reinserted as shown in figs. 22-23-24,


Fig.21


Fig.22


Fig.23


Fig.24

4.6 POSITIONING THE SPINAL BOARD INSIDE THE RESCUE BAG

4.6.1 With the patient immobilized on the spinal board

- Fully open the EVEREST bag by laying it on the ground and position the spinal board,
- fix the spinal board to the rescue bag by inserting the 6 buckles of the rescue bag into the slots of the spinal board, 4 of the upper “head-torso” section and 2 of the lower “legs” section (figs. 25, 26 and 27)


Fig.25


Fig.26


Fig.27

- close the zipper of the rescue bag starting from the feet of the patient (fig. 28),


Fig.28

d) close the rescue bag by overlapping the Velcro flaps (fig. 29) and fasten the external locking straps as shown in figs. 30, 31 and 32

Fig.29


Fig.30


Fig.31


Fig.32

e) place the mesh cover to repair the face of the injured person from the air and from any debris moved by the helicopter (fig. 33),


Fig.33

f) close the hood by pulling on the specific cords (fig. 34).


Fig.34

4.6.2 Without the patient

To reduce the preparation time of the kit during rescue operations, it is possible to prepare the spinal board, complete with head immobilizer, fixed in the rescue bag:

- a) open the spinal board without sliding the stiffening spars - paragraph 4.2 point a),
- b) if necessary, mount the head immobilizer on the spinal board as defined in the paragraph 4.3,
- c) place the spinal board in the sheet as defined in the paragraph 4.6.1 points a) and b) by inserting in the spinal board only the 4 buckles, placed in the upper "head-torso" section of the rescue bag (fig. 35),


Fig.35

- d) roll up the cover parts of the rescue bag, as shown in fig. 36, and fix them with the special velcro straps (figs. 37, 38 and 39),


Fig.36


Fig.37


Fig.38


Fig.39

e) fold the spinal board (fig. 40)

Do not leave materials inside (risk of damage to the spinal board and / or breakage of the rivets)


Fig.40

f) fasten the two external straps of the rescue bag (those of the lower part) as in fig. 41,


Fig.41

g) fold the rescue bag on itself, fasten the other two external straps (those on the upper part) and insert it into the backpack (fig. 42).


Fig.42

4.7 TRANSPORTING THE PATIENT

The patient, completely immobilized on the spinal board and inserted into the rescue bag, can be:

- manually transported by two operators using the fabric handles (fig. 43),
- attached to the helicopter winch,
- manually transported by two operators with the modular frame.


Fig.43

4.7.1 Connection with Patty hanging kit

- Connect the hanging kit to the fabric handles using the ten connectors and respecting the colors combination
- screw in the connector nuts completely (fig. 44)


Fig.44

CHAPTER 5

MAINTENANCE AND REPAIR

5.1 GENERAL INFORMATION

The MD is made of material which is highly resistant to wear and external agents. Despite this, the use conditions require maintenance operations and even repair in particular conditions.

Maintenance and repair operations must be recorded in a special register, an example of which is indicated in chapter 10.

5.2 MAINTENANCE

The maintenance operations that must be carried out by the user are, in order:

- a) cleaning: after every use, wash with warm (max. 40 °C) potable water, possibly adding a neutral detergent (e.g. Marseille soap). Rinse and leave to dry in a shaded place, away from direct sources of heat,
- b) disinfection, when deemed necessary or after prolonged inactivity (more than 21 days), preceded by cleaning: wipe with a clean cloth, soaked in a water solution of sodium hypochlorite (bleach) with a 0,1% concentration (1000 ppm). When blood or other bodily fluids are present, the recommended concentration of sodium hypochlorite is 0,5% (5000ppm).

5.3 REPAIR

Repair shall be strictly carried out by the manufacturer.
The user is allowed to replace the parts mentioned in paragraph 3.4.2 only, with new and genuine parts.

CHAPTER 6

6 STORAGE

After cleaning, any disinfection and drying operations, keep the product and its accessories in a place that is dry (relative humidity 40-90%), cool (temperature 5-40° C) and safe (avoid UV radiation), chemically neutral (absolutely avoid salty and/or acid environments), away from sharp edges, sources of heat, humidity, corrosive substances or other possible jeopardising conditions.

Do not store this device wet!

CHAPTER 7

CHECKS, INSPECTIONS AND SERVICING

7.1 CHECKS

Before and after every use, check the device and ensure that:

- it is appropriate for the use it was intended for,
- it has not been subjected to mechanical deformations and it neither has cracks nor shows signs of wear,
- all the rivets of the spinal board are intact,
- the textile parts do not have cuts, burns, chemical products residues, excessive hair, wear, in particular check areas at contact with the metal components,
- markings, including labels, are readable.

7.2 INSPECTIONS

The device must be inspected at least once a year, starting from the date of first use, by expert personnel trained and approved by KONG S.p.A. The date of the first use and the outcome of the inspection must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10.

7.3 SERVICING

The device must be serviced by KONG S.p.A., or by expert personnel trained and approved by KONG S.p.A. if:

- any malfunctions are detected,
- the result of the pre and post use check or of the inspections is negative.

The outcome of the servicing must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10.

The revised device is guaranteed for one year from the revision date.

CHAPTER 8

DURATION AND WARRANTY

8.1 DURATION

Carefully read the “Warnings and limitations of use” paragraph in chapter 2. The duration of this device depends on the outcome of the at least one-yearly inspection (paragraph 7.2) and possible servicing. The suspension, even if not used, has a maximum duration of 10 years from the date of manufacture indicated on the label.

Important: any devices that have not passed the pre-use, post-use and periodic inspection must be eliminated and made unusable.

8.2 DISPOSAL

Follow the rules in force in the country of use and the waste disposal procedures of the reference hospital facility for proper disposal.

8.3 WARRANTY

The manufacturer guarantees the compliance of the device with the regulations in force at the time of production. Warranty covering defects is limited to defects regarding raw materials and manufacturing defects. The warranty does not cover normal wear, oxidation, damage caused by non-compliant use and/or in competition, improper maintenance, transportation, preservation or storage. Warranty will be immediately deemed null and void should the device be modified or tampered with. Validity corresponds to the legal validity of the country where the device was sold, as from the date of sale, by the manufacturer. No claims shall be raised against the manufacturer once the aforementioned period has expired. Any repair or replacement request under warranty must be accompanied by a proof of purchase. Should the defect be recognised, the manufacturer undertakes to repair or, replace or refund the device at his discretion. The manufacturer shall not be held liable beyond the invoice price of the device in any case whatsoever.

8.4 LAW OBLIGATIONS

Professional and free time activities are often regulated by special national laws which can impose limits and/or obligations concerning the use of these devices. The user has the obligation to know and apply such laws which could provide for limits different from those indicated in this information.

CHAPTER 9

LABELLING AND SYMBOLS

9.1 MEDICAL DEVICE LABEL

This label features the KONG logo, a QR code, and various regulatory symbols. It includes the model name '855.00 - MOD. EVEREST', a UDI number '(01) 08023577044446(21) 1234561234', and the manufacturer's address: 'KONG S.p.A. Via XXV Aprile, 4 I - 23804 MONTE MARENZO (LC) - ITALY Tel. +39 0341 630506'. A factory icon is accompanied by the text 'ITALY 2022 / 01'.

This label features the KONG logo, a QR code, and various regulatory symbols. It includes the model name '855.11 - MOD. X-TRIM 2', a UDI number '(01) 08023577015057(21) 1234561234', and the manufacturer's address: 'KONG S.p.A. Via XXV Aprile, 4 I - 23804 MONTE MARENZO (LC) - ITALY Tel. +39 0341 630506'. A factory icon is accompanied by the text 'ITALY 2022 / 01'.

This label features the KONG logo, a QR code, and various regulatory symbols. It includes the model name '855.06 - MOD. HRP', a UDI number '(01) 08023577044606(21) 1234561234', and the manufacturer's address: 'KONG S.p.A. Via XXV Aprile, 4 I - 23804 MONTE MARENZO (LC) - ITALY Tel. +39 0341 630506'. A factory icon is accompanied by the text 'ITALY 2022 / 01'.

Any changes to the positioning of the symbols does not change the content thereof.

9.2 SEPARABLE COMPONENTS LABEL

Any changes to the positioning of the symbols does not change the content thereof.

9.3 SYMBOLS

Identification and traceability

	Medical Device
	Product identification code
	Unique Device Identifier
	Manufacturer's identification
	Unique serial number
	Country, year and month of production.

Data matrix UDI
 Readable with the Apps:
 iGepir - BarValid - Barcode Scanner

Information and safety

	Consult the user manual
	Compliant with Regulation (EU) 2017/745

CHAPTER 10

DOCUMENTS

10.1 INSPECTIONS AND SERVICING REGISTER

RESCUE BAG - Class I Medical Device				
REF.	855.00		MODEL	EVEREST
UDI-DI	08023577044446		SN	
Date of first use				
INSPECTION AND SERVICING				
DATE	I/S	Description	Outcome	Supervisor

SPINAL BOARD - Class I Medical Device

REF.	855.11	MODEL	X-TRIM 2
UDI-DI	08023577015057	SN	
Date of first use			

INSPECTION AND SERVICING

DATE	I/S	Description	Outcome	Supervisor

HEAD IMMOBILIZER – Class I Medical Device				
REF.	855.06	MODEL	HRP	
UDI-DI	08023577044606	SN		
Date of first use				
INSPECTION AND SERVICING				
DATE	I/S	Description	Outcome	Supervisor

10.2 MAINTENANCE AND REPAIR REGISTER

RESCUE BAG - Class I Medical Device				
REF.	855.00	MODEL	EVEREST	
UDI-DI	08023577044446	SN		
Date of first use				
MAINTENANCE AND REPAIR				
DATE	M/R	Description	Outcome	Supervisor

SPINAL BOARD – Class I Medical Device

REF.	855.11	MODEL	X-TRIM 2
UDI-DI	08023577015057	SN	
Date of first use			

MAINTENANCE AND REPAIR

DATE	M/R	Description	Outcome	Supervisor

HEAD IMMOBILIZER – Class I Medical Device			
REF.	855.06	MODEL	HRP
UDI-DI	08023577044606	SN	
Date of first use			

MAINTENANCE AND REPAIR				
DATE	M/R	Description	Outcome	Supervisor

10.3 EU DECLARATION OF CONFORMITY

The EU Declaration of Conformity is available for download at www.kong.it

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V08050199 – Stretchers for the transfer of patients - others

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Stretcher for transporting patients during rescue operations

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www.kong.it

11.1 APPLIED STANDARDS

- EN 597-1: 2015
Furniture - Assessment of the ignitability of mattresses and upholstered bed bases - Part 1: Ignition source smouldering cigarette
- EN 1865-1: 2010 + A1: 2015
Patient handling equipment used in road ambulances - Part 1: General stretcher systems and patient handling equipment
- EN 62366-1: 2015 + A1: 2020
Medical devices - Application of usability engineering to medical devices
- EN ISO 846:2019
Plastics - Evaluation of the action of microorganisms
- EN ISO 10993-1: 2020
Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
- EN ISO 10993-5: 2009
Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
- EN ISO 10993-23: 2021
Biological evaluation of medical devices - Part 23: Tests for irritation
- EN ISO 14971: 2019
Medical devices - Application of risk management to medical devices
- EN ISO 15223-1: 2021
Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements
- EN ISO 20417: 2021
Medical devices - Information to be supplied by the manufacturer.

11.2 STANDARDS USED AS REFERENCE

- EN ISO 13485: 2016
Medical devices - Quality management systems - Requirements for regulatory purposes
- EN 13718-1: 2014 + A1: 2020
Medical vehicles and their equipment - Air ambulances - Part 1: Requirements for medical devices used in air ambulances.

KONG s.p.a.

Via XXV Aprile, 4 - (zona industriale)

I - 23804 MONTE MARENZO (LC) - ITALY

Tel +39 0341630506 - Fax +39 0341641550

www.kong.it